ENG Bacteria-binding super-absorbent dressing

Device Description

Cutimed® Sorbion® Sorbact® is a sterile, bacteria and fungi binding, superabsorbent wound dressing. It consists of a green Cutimed® Sorbact® wound contact layer combined with a superabsorbent core and a white, water repellent backing, which prevents exudate strike through.

Cutimed® Sorbion® Sorbact® absorbs and retains exudate, thereby reducing the risk of maceration and enabling a moist wound environment. The high absorption capacity of Cutimed® Sorbion® Sorbact® helps decrease dressing change frequency. The dressing can be used in compression therapy.

Mode of Action

Cutimed® Sorbion® Sorbact® binds microorganisms, such as Staphylococcus aureus (including MRSA), Streptococcus species, Escherichia coli, Pseudomonas aeruginosa and Candida albicans, as shown in vitro. These microorganisms are removed from the wound each time the dressing is changed.

Intended Use

Cutimed® Sorbion® Sorbact® is intended for use in the management of clean, contaminated, colonized or infected wounds with high to excessive exudate levels, such as surgical wounds, traumatic wounds, pressure ulcers, diabetic ulcers and foot and leg ulcers.

Contraindications

There are no known contraindications to the use of Cutimed® Sorbion® Sorbact®. Attention is drawn to the warnings and precautions.

Warnings and Precautions

If the treated condition deteriorates, fails to improve or if a side effect is observed, consult a physician or an appropriate health care provider.

Do not use on eyes, mucous membranes or in wound cavities. The dressing swells considerably after fluid absorption.

Do not use on arterial bleeds or heavily bleeding wounds.

Do not use on patients with known hypersensitivity to any of the dressing components.

The green wound contact layer should always be applied in direct contact with the wound area to allow microorganisms to bind to the dressing.

Should the dressing dry or adhere to the wound, moisten or soak the dressing to assist removal and avoid disruption of the healing wound.

Do not use in combination with fatty products, such as ointments, creams and solutions, as they may decrease the binding of microorganisms.

Cutimed® Sorbion® Sorbact® is for single use only and should not be re-used. Re-use may lead to increased risk of infection or cross-contamination.

Do not re-sterilize.

Do not cut the dressing.

Do not use if peel pouch is already open or damaged as the sterility of the dressing is guaranteed only when the package is unopened and undamaged prior to use.

Prior to radiation therapy remove Cutimed® Sorbion® Sorbact®. A new dressing can be applied following treatment. Cutimed® Sorbion® Sorbact® is magnetic resonance imaging (MR) safe.

Instructions for Use

1. Prepare the wound according to local clinical practice. If desired, gently press a Cutimed® Sorbact® Round Swab against the moist wound surface to remove superficial microbial contamination

2. Select an appropriate dressing size for the wound. The wound pad should overlap the wound margins by at least 2 cm.

3. Remove the dressing from the pouch using an aseptic technique.

4. Apply the dressing. Ensure that the green side of the dressing comes into direct contact with the surface of the wound area. The white side of the dressing needs to face the caregiver.

5. Secure with tape at the dressing edges or cover with bandage. Ensure the dressing is able to expand under the fixation.

6. The dressing change frequency depends on exudate levels and the overall condition of the wound and surrounding skin. If saturated, the dressing shall be changed. Should the clinical condition allow, the dressing can be left in place for up to 7 days. The final decision on the therapy regimen needs to be taken by health care professionals.

Storage and Disposal

Keep dry and away from sunlight. Store below 40°C/104°F. Dispose of waste according to local environmental procedures.

Notice
Any serious incident that has occurred in relation to this device should be reported to BSN medical GmbH and the competent authority of your state.

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