Instructions for Use

Cutimed® Siltec® B

(ENG)

Product description

Cutimed® Siltec® B dressings are an assortment of silicone coated foam dressings. Cutimed® Siltec® B is a sterile, absorbent polyurethane foam dressing for traumatic dressing changes that contains super-absorbent strips which absorb and lock wound exudate. The wound contact surface is a perforated silicone layer that allows adherence to the peri-wound skin but not to the moist wound bed or to newly formed epithelial tissue. This minimises trauma and pain during dressing changes. The outer film is water-repellent, yet permeable to oxygen and vapour. The foam dressing and the super-absorbent strips are designed to absorb and to lock away excess wound fluid inside the dressing, promoting a moist wound environment and minimising the risk of maceration. Cutimed® Siltec® B is highly conformable and provides protection of the wound site. Cutimed® Siltec® B leaves no residues in the wound. The dressing reliably retains wound exudate, even under compression. If necessary the dressings can be cut to size. Cutimed® Siltec® B has an additional silicone adhesive border, that allows secure and gentle fixation of the dressing.

Indications

Cutimed® Siltec® B is indicated for exuding wounds with low to high exudate levels such as: venous and arterial ulcers, pressure ulcers, diabetic gangrenes, surgical incisions, skin grafts, lacerations or abrasions. Cutimed® Siltec® B is recommended for wounds with low to high exudate levels.

Contraindications

Do not use on patients with a known sensitivity to any of the product components.

Warnings and precautions

Cutimed® Siltec® B is packaged for single use. Do not re-use or re-sterilise as there is a risk of transmission of body fluids or contaminated tissue between patients. Do not use if pouch is open or damaged. Discard open or unused material. The wound should be inspected for signs of infection and treated according to clinical practice if required. For external use only. Do not use with oxidising solutions such as hypochlorite or hydrogen peroxide. Remove the dressing prior radiation therapy.

Directions for use

- Clean the wound according to good clinical practice.
- Make sure, the surrounding skin is dry.
- Choose an appropriate dressing size that covers the wound and overlaps the wound edges. **(A)**
- Open the pouch and take out the dressing. Ensure sterility.
- Remove protective film **(B)** and apply with the adherent side towards the wound. **(C)**
- Avoid stretching the dressing.
- Depending on wound location it may be appropriate to apply additional dressing fixation.
- If necessary, Cutimed® Siltec® B can be used in combination with other dressings e.g. Cutimed® Sorbact® (in case of infected wounds) or Cutimed® Gel (for dry, necrotic wounds).
- Cutimed® Siltec® B may be used on an infected wound in addition to antimicrobial therapy under medical supervision.
- Cutimed® Siltec® B can remain in place for up to seven days, depending on wound conditions and exudate levels. It is recommended to change the dressing every 24 hours initially, moving to less frequent changes as necessary.
- Carefully remove any additional dressing fixation and gently remove Cutimed® Siltec® B from the wound and discard.
- If required, cleanse the wound prior to application of a new dressing.

Storage

Cutimed® Siltec® B should be stored in dry conditions. Do not expose to direct sunlight.

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



The following artwork was approved:

Asset Name	JBR09149.01
Component Name	JBR09149.01
Project Title	C0797 - Siltec re-design IFUs
Revision	1
Reason for Revision	re-design

Asset Status

Date (GMT)	Signed By	Role	Action	Comment
2016-06-16 16:56:10.0	Anja Spruth	Artwork Management	Artwork Submitted for Approval	I am the author of this document
2016-07-01 12:52:29.792	Helgo Schieber	Packaging Development	Artwork Approved : Electronic Signature	I attest to the accuracy and integrity of this document
2016-07-01 12:29:39.633	Nicole Bartels	Global Commercial	Artwork Approved : Electronic Signature	I attest to the accuracy and integrity of this document
2016-07-01 11:55:51.875	Janina Hinck	Regulatory Affairs HAM	Artwork Approved : Electronic Signature	I attest to the accuracy and integrity of this document
2016-07-01 13:05:08.0	Anja Spruth	Artwork Management	Final Sign Off	I attest to the accuracy and integrity of this document

THIS REPORT IS CREATED FROM ELECTRONIC SIGNATURE ACTIVITY WITH KALLIK AMS 360

CREATED ON 2016-07-01 13:05:08.0

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BSN medical GmbH