



medi

circaid[®] juxtalite[®]

Adjustable compression device for treating venous diseases.

medi. I feel better.

circaid[®] juxtalite[®]

Adjustable compression device for treating venous diseases.

medi Tipp*

Due to product features we recommend circaid juxtalite for:

- Women / men
- Patients who have challenges with medical compression stockings
- 24-hour compression therapy
- Mobile / immobile patients

Indications

- Chronic venous insufficiency
- Varicose veins
- Post Thrombotic Syndrome
- Venous Stasis Ulcer
- Post Sclerotherapy
- Dependent oedema

Product features

- Easy handling – simple donning and doffing
- circaid Built-In-Pressure System (BPS): adjustable, controllable and re-adjustable compression
- Breathe-O-Prene: breathable material with high working pressure
- Latex-free, odour-inhibiting material with anti-bacterial effect
- Re-adjustability of the individual bands avoids pressure loss and thus sliding down of the garment^{1,2}
- **New:** Also available with new compressive undersock
 - Applies graduated compression to the foot and ankle
 - Graduated compression continues into the leg portion of the undersock

Benefits

- Offering an alternative therapy concept for patients who are not able to apply or tolerate ulcer compression stockings
- Faster healing rates of venous leg ulcer^{1,2}
- Promotes self-management thanks to easy handling¹



circaid juxtalite with circaid compression ankle



circaid juxtalite with circaid compressive undersock

Note


For more information and handling videos please visit: www.medi.biz/jl

* Please consider, that this recommendation does not replace the diagnosis and consultation of a medical healthcare professional.

Studies

¹ Mosti G et al. Adjustable compression wrap devices are cheaper and more effective than inelastic Bandages for venous leg ulcer healing. A Multicentric Italian Randomized Clinical Experience. *Phlebology* 2020;35(2):124-133.

² Stather PW et al. Review of adjustable velcro wrap devices for venous ulceration. *In Wound J* 2019;16(4):903-908.



Circumferences in cm								
	S	M	M-X	L	L-X	XL	XL-X	XXL
c C	26-36	33-44	44-54	41-51	51-61	37-48	48-61	53-64
c B	19-22	22-27	22-27	25-32	25-32	30-37	30-37	35-42

Actual garment length in cm		
	Short	Long
∠ B-D	28	33
Sizing measurement in cm		
∠ floor - knee crease	< 44	> 44

Size / Article No.	circaid juxtalite		circaid juxtalite with circaid compressive undersock	
	Short	Long	Short	Long
S	JU15702000	JU15002000	JU1M702000	JU1M002000
M	JU15703000	JU15003000	JU1M703000	JU1M003000
M-X	JU15713000	JU15013000	JU1M713000	JU1M013000
L	JU15704000	JU15004000	JU1M704000	JU1M004000
L-X	JU15714000	JU15014000	JU1M714000	JU1M014000
XL	JU15705000	JU15005000	JU1M705000	JU1M005000
XL-X	JU15715000	JU15015000	JU1M715000	JU1M015000
XXL	JU15706000	JU15006000	JU1M706000	JU1M006000

Your order includes

- 1 inelastic adjustable compression garment
- 1 circaid Built-In-Pressure System (BPS) guide card
- 1 pair of circaid compression anklets
- 1 circaid undersleeve lower leg


Your order includes

- 1 inelastic adjustable compression garment
- 1 circaid Built-In-Pressure System (BPS) guide card
- 1 pair of circaid compressive undersocks moderate

circaid® compressive undersock

Easy-to-apply liner with targeted compression for foot and ankle.

- Applies graduated compression to the foot and ankle
- Graduated compression continues into the leg portion of the undersock
- Soft material for pleasant wearing comfort
- Available in 15 – 25 mmHg (moderate compression) and 25 – 35 mmHg (firm compression)
- Available in standard and large



Circumferences in cm		
	Standard	Large
cA1	19-32	24-41
cB	19-27	25-42
cC	26-54	41-64

Size	Strength	Article No.
Standard	Moderate	JU5M001 000
Large	Moderate	JU5M002 000
Standard	Firm	JU5F001 000
Large	Firm	JU5F002 000

