

# TECHNICAL DATA SHEET



## Description

### Reference

#### ENFit® feeding connectors

50-6112	CORFLO* PEG Replacement Feeding Adapter Fr 12 (5x1)
50-6116	CORFLO* PEG Replacement Feeding Adapter Fr 16 (5x1)
50-6120	CORFLO* PEG Replacement Feeding Adapter Fr 20 (5x1)

#### NON-ENFit® feeding connectors

30-6112	CORFLO* PEG Replacement Feeding Adapter Fr 12 (5x1)
30-6116	CORFLO* PEG Replacement Feeding Adapter Fr 16 (5x1)
30-6120	CORFLO* PEG Replacement Feeding Adapter Fr 20 (5x1)

This CORFLO PEG feeding connector is used to replace a broken or used CORFLO\* PEG feeding connector. The size of the feeding connector should match the size of the CORFLO\* PEG tubing. The products with item#50-61xx have a male ENFit® connector as feeding port.

The products with item# 30-61xx have a female Luer Lock connection

### Kit composition:

- 1 Feeding connector
- 1 Gauze 10x10cm

## Indication

This allows to replace the broken or damaged component of the PEG and prevents the patient from undergoing an endoscopic intervention for replacement of the complete device.

The replacement of the feeding connector can, when the caregiver is properly trained and instructed, be done in the alternate care segment.

## Counter indication

Do not replace the feeding connector when the appropriate size of feeding connector is not available or when the PEG tubing itself is damaged.

## Main materials

Feeding connector: Copolyester

Closure caps and strap: Copolyester/Thermoplastic Elastomer  
Not formulated with Natural Rubber Latex.

Not formulated with DEHP / DBP / BBP / Bisphenol A.



## Packaging

Each feeding connectors is packed together with the gauze in a plastic pouch.

5 plastic pouches are paked in a carton box. a lamps are packed in a plastic box.

Size: (LxWxD) 29,2 x 13,0 x 7,9 cm

Box weight: ... g

Bar coding: GS1-128 symbology, linear, on unit packaging and carton box.

Item #	Single Unit	Carton box
50-6112	(01) 10680651462128	(01) 30680651462122
50-6116	(01) 10680651462135	(01) 30680651462139
50-6120	(01) 10680651462142	(01) 30680651462146
30-6112	(01) 10680651458916	(01) 30680651458910
30-6116	(01) 10680651458923	(01) 30680651458927
30-6120	(01) 10680651458930	(01) 30680651458934

## Manufacturing

Manufactured Avent Tucson.

## Regulatory information

Product CE marked as per 93/42/EEC Directive.

Class of the device: I.

Notified body: BSI (0086)

## Sterilisation

These devices are non-sterile

## Storage

Store in a dry and cool place, away from sources of heat and radiation. Keep as much as practicably possible in its carton box.

## Shelf life

5 years, from date of manufacturing.

\*Registered trademark or trademark of Halyard Health Inc.

The dimensions and properties listed above can vary within pre-established specifications. This document was created using the most recent information. In the interest of continuous improvement, the characteristics of the product may change without prior notice.