


Standard Single Use Syringes



Single Use Syringes Luer Lock \geq 10ml



1. The narrow silicone seal minimizes the contact area between the silicone and the solution thereby reducing solid particle contamination. The silicone seal provides easy and accurate dosage management.
2. Luer-Lock Nozzle The dual scale provides easy and precise reading irrespective of the handling of the syringe. Plunger stop prevents accidental removal of the plunger.
3. The plunger requires despite of reduced amount of silicone oil a minimum use of force due to smooth plunger action. Finger grip of an ergonomic design makes the syringe easy to handle.

REF	Product Description	Packaging Units
62.6711	10 ML CODAN Luer-Lock 0,2 grad Bulk Non-Sterile	1800/28800
62.6716	10 ML CODAN Luer-Lock 0,2 grad Sterile Bag	400//800/12800
62.6706	10 ML CODAN Luer-Lock 0,2 grad	100/1200/16200
62.6703	10 ML CODAN Luer-Lock 0,5 grad	100/1200/16200
62.6708	10 ML CODAN Luer-Lock 1,0 grad	100/1200/16200
62.7709	20 ML CODAN Luer-Lock 1,0 grad Bulk Non-Sterile	900/14400
62.7716	20 ML CODAN Luer-Lock 0,5 grad Sterile Bag	280/560/8960
62.7711	20 ML CODAN Luer-Lock 0,5 grad	50/800/12800
62.7704	20 ML CODAN Luer-Lock 1,0 grad	50/800/12800
62.7705	20 ML CODAN Luer-Lock 2,0 grad	50/800/12800